



May 2011

Chlamydia trachomatis and Neisseria gonorrhoeae - Specimen Collection Guidelines

To Health Care Providers:

Update:

- Vaginal swabs are not accepted for CT/GC NAA testing at the OAHPP – Public Health Laboratories.
• Testing for Neisseria gonorrhoeae (GC) as a single test by Nucleic Acid Amplification is no longer performed.

Table with 3 columns: Collection Kit, Specimen, When to use...... Rows include Unisex Swab Specimen Collection Kit and Urine Specimen Collection Kit.

Table with 3 columns: Collection Kit, Specimen, When to use...... Rows include Chlamydia trachomatis and Neisseria gonorrhoeae. Includes a text block at the bottom regarding Test of Cure for CT/GC.

For Further Information:

- OAHPP Customer Service Centre toll free at 1-877-604-4567 or 416-235-6556.
• For the OAHPP Specimen Collection Guide and previous Labstracts, refer to http://www.oahpp.ca/publichealthlaboratories.php